

# IoT Based Automated Myo Neuro Stimulator For Neuro Paralytic Patients

Submitted By:

**Team : VARAHI**

**SRI SAIRAM INFOTECH**

**Mentor:** Vice President - Sales & Marketing  
AKAS Medical, Chennai



Mr.R.Dhanasekaran



Mr.M.ArunPrakash



Mr.R.Hemanath



Mr.K.K.ArunRitik

Ms.K.Vidhya ,  
Professor,KPRIET



Dr.P.Rathan Pandian  
Community Physician

# Objective of Innovation

- **Bring back Near-Normal Life(Rehabilitation) of Neuro -Paralytic Patients**
- **Self motivated, Regularized exercises**
- **No manual assistance-Normal Routine Activities**

# Problem Addressed

We all know someone - a brother, sister, friend, neighbor, or colleague -**living with paralysis**. They **aren't strangers**, But

- ✓ **They are mostly unnoticed and separated**
- ✓ **Need maximum care with compulsory manual assistance for rehabilitation**


# Motivation

## Paralytic Patient Care – Faster and Self Motivated Recovery of Paralytic Patients

Paralysis Prevalence and Health Disparities Survey (PPHDS)-2020:

- ✓ Approximately 2.5% of total population are living with paralysis
- ✓ 1.8 % of paralytic patients are hopeless and completely dependent

# Neuro Pralysis



- ✓ Loss of strength in and control over a muscle or group of muscles .
- ✓ More likely due to **damage of chain of nerve cells** that runs from the body part to brain and back again.

(Nerve cells deliver the signals for our muscles to move)

# Neuro Paralysis

## Etiological Cause of Neuro Paralysis

- ✓ Guillain-Barré syndrome
- ✓ Stroke
- ✓ Myasthenia Gravis
- ✓ Botulism
- ✓ Drugs
- ✓ Head Injury/Trauma
- ✓ Infections

## Clinical Features of Neuro Paralysis

- ✓ Muscle Weakness
- ✓ Stiffness of one/both sides of body
- ✓ Muscle Spasticity/Permanently Contracted Muscle
- ✓ Poor Fine Motor Skills
- ✓ Trouble Walking
- ✓ Poor Balance
- ✓ Cranial/Spinal Nerve dysfunction

# Patient Rehabilitation

## Exercise -Foundation of the healing and recovery process

- ✓ **Active rehab** exercise -Movement on our own
- ✓ **Passive exercise** requires **assisting movement**
  - ✓ **Require help of a caregiver, therapist, or family member** to move the affected arm or leg
- Needs pain-free ,comfortable and friendly exercise platform
- Needs more positive and energized enforcement (Self Interest)



# EXISTING SOLUTION



# 1.Exercise-Manually Assisted




## 2. Electrical Stimulation for Stroke Rehabilitation



➤ It is used for **paralyzed lower limbs after stroke** in preventing forefoot dropping

### 3. Belt Electrode Skeletal Muscle Electrical Stimulation (B-SES)



➤ B-SES - stimulate all the muscles in the lower extremities



## 4.Functional Electrical Stimulation(FES) Devices



A. NESS L300 and  
wireless foot switch  
for gait

B. NESS H200, which is  
worn over the  
paralyzed arm and  
hand for upper limb  
motion

# Proposed Innovation


## “IoT Based Myo Neuro Stimulator for Neuro Paralytic Patients”

### Goal:


To **recover** motor loss and sensory loss in the affected parts of paralytic patients through exercise,

- ✓ In absence of care givers
- ✓ Patient Convenient Time

# BLOCK DIAGRAM




# Work Flow





# SCHEMATIC DIAGRAM OF INNOVATION



# COMPONENETS

- A Microcontroller (1)
- A Real Time Controller (2)
- A DC motor and Driver module (3)
- A Li-Po battery (4)
- A SD Card Module (5)
- A Speaker (6)
- A Voice Recognition Module (7)
- A GSM module (8) and
- A Blink Sensor (9)



# Novelty of Proposed Innovation



## Facilitates:

### 1. Safety Features – Limit Switch Controlled

- ✓ Controlled Angular Movement
- ✓ Position Based and Timer based – Auto Off

### 2. Regularized - Timer Based Patient Alert

### 3. Interest based Activation

### 4. Patient Convenient Control

- ✓ Voice Control
- ✓ Remote Control



# Novelty of Proposed Innovation

## 5. Selected Activation

Convenient based body part activation(Based on pain and cramp)

## 6. Relaxed Journey

Mentally Relaxed, Motivated Journey –Musical System Enabled

## 7. Emergency Communication

To care givers- GSM Supported

## 8. Intelligent Evaluation and Display

Display the count and period of exercise

# Benefits of Proposed Innovation

- Prevents patients from atrophy or muscle loss due to physical inactivity
- Facilitates - Regularized exercises ,even in the absence of their care takers
- Patients interest based body part activation
- Patient convenient

# Benefits of Proposed Innovation

- Eliminates the need of long time assistance
- Facilitates adjustable arm and leg holder
- Automated process of exercise based on patients interest
- Remote and voice recognition system
- **GSM Support** -To contact for emergency needs
- **Music System** -Maintain patients' attention and interest until completion



# EXISTING Vs PROPOSED

# Existing Vs Proposed

- **Non-Electrode System:** Non placement of electrodes on the human body
  - No direct exposure to electric voltage across the body
- **Fully Automated:** No Human presence (Care taker) is warranted always
- **Patients Convenient and interest:** Start and Stop based on interest
- **Mode of Activation/Control :** **Complete controlled automation** through **Remote or Voice**

# EXISTING Vs PROPOSED

- **Aligned with the guidelines** of **Physiotherapist**
- **SMS Communication**  
Single button based or voice controlled SMS alert system to care taker of respective patients.
- **Psychological encouragement or stimulation**  
Facilitated through **music system as part of this invention.**
- **Cost** :The invention would be cost effective and affordable to lower class and middle class people.

# FUTURE EXTENSION

- Communication support for fully paralytic patients through Eye Blink Sensor Module.
- Visualization of Periodic Improvement

# MILESTONES REACHED

# 1. Patent

## Filed and Published patent for the proposed invention



Controller General of Patents, Designs and Trademarks  
Department of Industrial Policy and Promotion  
Ministry of Commerce and Industry

Application Details	
APPLICATION NUMBER	202041031299
APPLICATION TYPE	ORDINARY APPLICATION
DATE OF FILING	22/07/2020
APPLICANT NAME	1 . K. Vidhya 2 . Dr. P. Rathan 3 . M. Premkumar Balaji D
TITLE OF INVENTION	IOT BASED AUTOMATED MYONEURO STIMULATOR FOR NEURO PARALYTIC PATIENTS
FIELD OF INVENTION	ELECTRONICS
E-MAIL (As Per Record)	
ADDITIONAL-EMAIL (As Per Record)	vidhya.k@kpriet.ac.in
E-MAIL (UPDATED Online)	
PRIORITY DATE	
REQUEST FOR EXAMINATION DATE	22/07/2020
PUBLICATION DATE (U/S 11A)	07/08/2020

### Application Status

# 2.Startup Registration(MSME)

## 2.Startup Registered under MSME

		<p>भारत सरकार Government of India सूक्ष्म, लघु एवं मध्यम उद्यम मंत्रालय Ministry of Micro, Small and Medium Enterprises</p>																							
<p><b>UDYAM REGISTRATION CERTIFICATE</b></p>																									
 <p>Our small hands to make you LARGE</p>																									
TYPE OF ENTERPRISE	MICRO		MANUFACTURING																						
UDYAM REGISTRATION NUMBER	UDYAM-TN-03-0001862																								
NAME OF ENTERPRISE	SRI SAIRAM INFOTECH																								
OFFICAL ADDRESS OF ENTERPRISE	<table><tr><td>Flat/Door/Block No.</td><td>2ND FLOOR</td><td>Name of Premises/ Building</td><td>-</td></tr><tr><td>Village/Town</td><td>SIDDHAPUDUR</td><td>Block</td><td>SIDDHAPUDUR</td></tr><tr><td>Road/Street/Lane</td><td>52-53 NGN STREET</td><td>City</td><td>SIDDHAPUDUR</td></tr><tr><td>State</td><td>TAMIL NADU</td><td>District</td><td>COIMBATORE , Pin 641044</td></tr><tr><td>Mobile</td><td>9865511224</td><td>Email:</td><td>vidhyasenthilkumar1@gmail.com</td></tr></table>					Flat/Door/Block No.	2ND FLOOR	Name of Premises/ Building	-	Village/Town	SIDDHAPUDUR	Block	SIDDHAPUDUR	Road/Street/Lane	52-53 NGN STREET	City	SIDDHAPUDUR	State	TAMIL NADU	District	COIMBATORE , Pin 641044	Mobile	9865511224	Email:	vidhyasenthilkumar1@gmail.com
Flat/Door/Block No.	2ND FLOOR	Name of Premises/ Building	-																						
Village/Town	SIDDHAPUDUR	Block	SIDDHAPUDUR																						
Road/Street/Lane	52-53 NGN STREET	City	SIDDHAPUDUR																						
State	TAMIL NADU	District	COIMBATORE , Pin 641044																						
Mobile	9865511224	Email:	vidhyasenthilkumar1@gmail.com																						
DATE OF INCORPORATION / REGISTRATION OF ENTERPRISE	13/08/2020																								



# 3.Prototype

## 3.Completion of Prototype Development(80%)



# MILESTONES REACHED



# QUERIES